



PLXNB1 Antibody, FITC conjugated

Product Code	CSB-PA018222LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43157
Immunogen	Recombinant Human Plexin-B1 protein (664-837AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plexin-B1 (Semaphorin receptor SEP), PLXNB1, KIAA0407 PLXN5 SEP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PLXNB1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.