



SELP Antibody, Biotin conjugated

Product Code	CSB-PA020978LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16109
Immunogen	Recombinant Human P-selectin protein (543-767AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	P-selectin (CD62 antigen-like family member P) (Granule membrane protein 140) (GMP-140) (Leukocyte-endothelial cell adhesion molecule 3) (LECAM3) (Platelet activation dependent granule-external membrane protein) (PADGEM) (CD antigen CD62P), SELP, GMRP GRMP
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	SELP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.