



HEXA Antibody, Biotin conjugated

Product Code	CSB-PA12259D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06865
Immunogen	Recombinant Human Beta-hexosaminidase subunit alpha protein (89-529AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-hexosaminidase subunit alpha (EC 3.2.1.52) (Beta-N-acetylhexosaminidase subunit alpha) (Hexosaminidase subunit A) (N-acetyl-beta-glucosaminidase subunit alpha), HEXA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	HEXA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.