



DOK7 Antibody, FITC conjugated

Product Code	CSB-PA623281LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q18PE1
Immunogen	Recombinant Human Protein Dok-7 protein (416-504AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Dok-7 (Downstream of tyrosine kinase 7), DOK7, C4orf25
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DOK7