



GPR158 Antibody, FITC conjugated

Product Code	CSB-PA735975LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T848
Immunogen	Recombinant Human Probable G-protein coupled receptor 158 protein (882-1003AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable G-protein coupled receptor 158, GPR158, KIAA1136
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPR158
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.