



SENP7 Antibody, HRP conjugated

Product Code	CSB-PA805340LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQF6
Immunogen	Recombinant Human Sentrin-specific protease 7 protein (128-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sentrin-specific protease 7 (EC 3.4.22.-) (SUMO-1-specific protease 2) (Sentrin/SUMO-specific protease SENP7), SENP7, KIAA1707 SSP2 SUSP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SENP7