



SENP7 Antibody, Biotin conjugated

Product Code	CSB-PA805340LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQF6
Immunogen	Recombinant Human Sentrin-specific protease 7 protein (128-366AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sentrin-specific protease 7 (EC 3.4.22.-) (SUMO-1-specific protease 2) (Sentrin/SUMO-specific protease SENP7), SENP7, KIAA1707 SSP2 SUSP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SENP7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.