



IL17F Antibody, Biotin conjugated

Product Code	CSB-PA842733LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PD4
Immunogen	Recombinant Human Interleukin-17F protein (31-163AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-17F (IL-17F) (Cytokine ML-1), IL17F
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IL17F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.