



OTUD5 Antibody, Biotin conjugated

Product Code	CSB-PA850298LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96G74
Immunogen	Recombinant Human OTU domain-containing protein 5 protein (51-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	OTU domain-containing protein 5 (EC 3.4.19.12) (Deubiquitinating enzyme A) (DUBA), OTUD5
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OTUD5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.