



MYO5A Antibody, HRP conjugated

Product Code	CSB-PA857430LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4I1
Immunogen	Recombinant Human Unconventional myosin-Va protein (1087-1220AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Unconventional myosin-Va (Dilute myosin heavy chain, non-muscle) (Myosin heavy chain 12) (Myosin-12) (Myoxin), MYO5A, MYH12
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MYO5A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.