



EIF1 Antibody, HRP conjugated

Product Code	CSB-PA007503LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41567
Immunogen	Recombinant Human Eukaryotic translation initiation factor 1 protein (2-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 1 (eIF1) (A121) (Protein translation factor SUI1 homolog) (Sui1iso1), EIF1, SUI1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EIF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.