



EIF1 Antibody, FITC conjugated

Product Code	CSB-PA007503LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41567
Immunogen	Recombinant Human Eukaryotic translation initiation factor 1 protein (2-113AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Eukaryotic translation initiation factor 1 (eIF1) (A121) (Protein translation factor SUI1 homolog) (Sui1iso1), EIF1, SUI1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	EIF1