



RAMP2 Antibody, FITC conjugated

Product Code	CSB-PA019305LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60895
Immunogen	Recombinant Human Receptor activity-modifying protein 2 protein (43-145AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor activity-modifying protein 2 (Calcitonin-receptor-like receptor activity-modifying protein 2) (CRLR activity-modifying protein 2), RAMP2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RAMP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.