



RASAL1 Antibody, Biotin conjugated

Product Code	CSB-PA019351LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95294
Immunogen	Recombinant Human RasGAP-activating-like protein 1 protein (698-804AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RasGAP-activating-like protein 1 (RAS protein activator like 1) (Ras GTPase-activating-like protein), RASAL1, RASAL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RASAL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.