



TOM1L1 Antibody, HRP conjugated

Product Code	CSB-PA024044LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75674
Immunogen	Recombinant Human TOM1-like protein 1 protein (296-476AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TOM1-like protein 1 (Src-activating and signaling molecule protein) (Target of Myb-like protein 1), TOM1L1, SRCASM
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TOM1L1