



ALG11 Antibody, HRP conjugated

Product Code	CSB-PA647192LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2TAA5
Immunogen	Recombinant Human GDP-Man:Man(3)GlcNAc(2)-PP-Dol alpha-1,2-mannosyltransferase protein (416-480AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDP-Man:Man(3)GlcNAc(2)-PP-Dol alpha-1,2-mannosyltransferase (EC 2.4.1.131) (Asparagine-linked glycosylation protein 11 homolog) (Glycolipid 2-alpha-mannosyltransferase), ALG11, GT8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ALG11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.