



SLC25A24 Antibody, HRP conjugated

Product Code	CSB-PA743993LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NUK1
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-1 protein (1-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein SCaMC-1 (Mitochondrial ATP-Mg/Pi carrier protein 1) (Mitochondrial Ca(2+)-dependent solute carrier protein 1) (Small calcium-binding mitochondrial carrier protein 1) (Solute carrier family 25 member 24), SLC25A24, APC1 MCSC1 SCAMC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A24