



SLC25A24 Antibody, HRP conjugated

Product Code	CSB-PA743993LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6NUK1
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-1 protein (1-197AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein SCaMC-1 (Mitochondrial ATP-Mg/Pi carrier protein 1) (Mitochondrial Ca(2+)-dependent solute carrier protein 1) (Small calcium-binding mitochondrial carrier protein 1) (Solute carrier family 25 member 24), SLC25A24, APC1 MCSC1 SCAMC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A24
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.