



RXFP2 Antibody, Biotin conjugated

Product Code	CSB-PA819908LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXD0
Immunogen	Recombinant Human Relaxin receptor 2 protein (33-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Relaxin receptor 2 (G-protein coupled receptor 106) (G-protein coupled receptor affecting testicular descent) (Leucine-rich repeat-containing G-protein coupled receptor 8) (Relaxin family peptide receptor 2), RXFP2, GPR106 GREAT LGR8
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RXFP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.