



RAB3GAP2 Antibody, HRP conjugated

Product Code	CSB-PA872480LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2M9
Immunogen	Recombinant Human Rab3 GTPase-activating protein non-catalytic subunit protein (277-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab3 GTPase-activating protein non-catalytic subunit (RGAP-iso) (Rab3 GTPase-activating protein 150 kDa subunit) (Rab3-GAP p150) (Rab3-GAP150) (Rab3-GAP regulatory subunit), RAB3GAP2, KIAA0839
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RAB3GAP2