




Acetyl-HIST1H1A (K21) Antibody

Product Code	CSB-PA010376PA21acHU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02539
Immunogen	Peptide sequence around site of Acetyl-Lys (21) derived from Human Histone H1.1
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H1.1 (Histone H1a), HIST1H1A, H1F1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIST1H1A

Image



Immunofluorescent analysis of HepG2 cells (sodium butyrate, 30 mM, 4h) using CSB-PA010376PA21acHU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)