





HIST1H1B (Ab-17) Antibody

Product Code	CSB-PA010377PA17nphHU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16401
Immunogen	Peptide sequence around site of Ser (17) derived from Human Histone H1.5
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histone H1.5 (Histone H1a) (Histone H1b) (Histone H1s-3), HIST1H1B, H1F5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HIST1H1B

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA010377PA17nphHU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA010377PA17nphHU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.