




HIST1H1B (Ab-66) Antibody

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA010377PA66nforHU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P16401 |
| Immunogen | Peptide sequence around site of Lys (66) derived from Human Histone H1.5 |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IHC; Recommended dilution: IHC:1:20-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Histone H1.5 (Histone H1a) (Histone H1b) (Histone H1s-3), HIST1H1B, H1F5 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | HIST1H1B |

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA010377PA66nforHU at dilution of 1:100