




AFF2 Antibody

Product Code	CSB-PA001413LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51816
Immunogen	Recombinant Human AF4/FMR2 family member 2 protein (575-705AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AF4/FMR2 family member 2 (Fragile X E mental retardation syndrome protein) (Fragile X mental retardation 2 protein) (FMR2P) (Protein FMR-2) (Protein Ox19), AFF2, FMR2 OX19
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AFF2

Image



Western Blot

Positive WB detected in: A549 whole cell lysate

All lanes: AFF2 antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 145, 141, 48, 144, 106 kDa

Observed band size: 145 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.