






Recombinant Mouse Serine protease inhibitor A3N (Serpina3n)

Product Code	CSB-EP835578MO(M)
Abbreviation	Recombinant Mouse Serpina3n protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91WP6
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SFPDGTLGMDAAVQEDHDNGTQLDSLTLASINTDFAFSLYKELVLKNPDKNIVF SPLSISAALAVMSLGAKGNTLEEILEGLKFNLTTETSEADIHQGFGHLLQRLNQPK DQVQISTGSALFIEKRQQILTEFQEKAFTLYQAEAFADTFQPRQAKKLINDYV RKQTQGMKELVSDLDKRTLMLVNYIYFKAKWKVPFDPLDTFKSEFYAGKRR PVIVPMMSMEDLTPYFRDEELSCVTVELKYTGNASALFILPDQGRMQQVEAS LQPETLRKWKNSLKPRMIDELHLPKFISTDYSLEDVLSKLGIREVFSTQADLSAI TGTKDLRVSQVVHKAFLDVAETGTEAAAATGVKFPMSAKLYPLTVYFNRPFLI MIFDTETEIAPIAKIANPK
Source	E.coli
Target Names	Serpina3n
Protein Names	Recommended name: Serine protease inhibitor A3N Short name= Serpin A3N
Expression Region	21-418aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	60.8kDa
Protein Length	full length protein(mutation)
Image	



We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.