



CD20 Recombinant Monoclonal Antibody

Product Code	CSB-RA015007A0HU
Abbreviation	B-lymphocyte antigen CD20
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11836
Immunogen	A synthesized peptide
Species Reactivity	Human
Tested Applications	ELISA
Relevance	This protein may be involved in the regulation of B-cell activation and proliferation.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Rabbit IgG in 10mM phosphate buffered saline , pH 7.4, 150mM sodium chloride, 0.05% BSA, 0.02% sodium azide and 50% glycerol.
Purification Method	Affinity-chromatography
Isotype	Rabbit IgG
Clonality	Monoclonal
Alias	B-lymphocyte antigen CD20, B-lymphocyte surface antigen B1, Bp35, Leukocyte surface antigen Leu-16, Membrane-spanning 4-domains subfamily A member 1, CD20, MS4A1, CD20
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MS4A1
Clone No.	2C6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.