



AGRN Antibody, HRP conjugated

Product Code	CSB-PA001461LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00468
Immunogen	Recombinant Human Agrin protein (968-1130AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Agrin [Cleaved into: Agrin N-terminal 110 kDa subunit; Agrin C-terminal 110 kDa subunit; Agrin C-terminal 90 kDa fragment (C90); Agrin C-terminal 22 kDa fragment (C22)], AGRN, AGRIN
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	AGRN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.