



CCS Antibody, Biotin conjugated

Product Code	CSB-PA004854LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14618
Immunogen	Recombinant Human Copper chaperone for superoxide dismutase protein (145-274AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Copper chaperone for superoxide dismutase (Superoxide dismutase copper chaperone), CCS
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CCS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.