



CLN5 Antibody, HRP conjugated

Product Code	CSB-PA005567LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75503
Immunogen	Recombinant Human Ceroid-lipofuscinosis neuronal protein 5 protein (42-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ceroid-lipofuscinosis neuronal protein 5 (Protein CLN5) [Cleaved into: Ceroid-lipofuscinosis neuronal protein 5, secreted form], CLN5
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CLN5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.