



CRYAA Antibody, Biotin conjugated

Product Code	CSB-PA006007LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02489
Immunogen	Recombinant Human Alpha-crystallin A chain protein (98-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-crystallin A chain (Heat shock protein beta-4) (HspB4) [Cleaved into: Alpha-crystallin A(1-172); Alpha-crystallin A(1-168); Alpha-crystallin A(1-162)], CRYAA, CRYA1 HSPB4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CRYAA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.