



CYBA Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA006323LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P13498 |
| Immunogen | Recombinant Human Cytochrome b-245 light chain protein (128-195AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Cytochrome b-245 light chain (Cytochrome b(558) alpha chain) (Cytochrome b558 subunit alpha) (Neutrophil cytochrome b 22 kDa polypeptide) (Superoxide-generating NADPH oxidase light chain subunit) (p22 phagocyte B-cytochrome) (p22-phox) (p22phox), CYBA |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | CYBA |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |