



CYBA Antibody, FITC conjugated

Product Code	CSB-PA006323LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13498
Immunogen	Recombinant Human Cytochrome b-245 light chain protein (128-195AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome b-245 light chain (Cytochrome b(558) alpha chain) (Cytochrome b558 subunit alpha) (Neutrophil cytochrome b 22 kDa polypeptide) (Superoxide-generating NADPH oxidase light chain subunit) (p22 phagocyte B-cytochrome) (p22-phox) (p22phox), CYBA
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CYBA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.