




DHX15 Antibody

Product Code	CSB-PA006870LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43143
Immunogen	Recombinant Human Pre-mRNA-splicing factor ATP-dependent RNA helicase DHX15 protein (371-502AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-mRNA-splicing factor ATP-dependent RNA helicase DHX15 (EC 3.6.4.13) (ATP-dependent RNA helicase #46) (DEAH box protein 15), DHX15, DBP1 DDX15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DHX15

Image



Western Blot

Positive WB detected in: 293 whole cell lysate

All lanes: DHX15 antibody at 3.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 91 kDa

Observed band size: 91 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.