



ESYT2 Antibody, Biotin conjugated

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA008279LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A0FGR8 |
| Immunogen | Recombinant Human Extended synaptotagmin-2 protein (567-736AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Extended synaptotagmin-2 (E-Syt2) (Chr2Syt), ESYT2, FAM62B KIAA1228 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ESYT2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |