



FFAR1 Antibody, FITC conjugated

Product Code	CSB-PA008604LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14842
Immunogen	Recombinant Human Free fatty acid receptor 1 protein (143-178AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Free fatty acid receptor 1 (G-protein coupled receptor 40), FFAR1, GPR40
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FFAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.