



FGFRL1 Antibody, HRP conjugated

Product Code	CSB-PA008650LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N441
Immunogen	Recombinant Human Fibroblast growth factor receptor-like 1 protein (403-470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor receptor-like 1 (FGF receptor-like protein 1) (FGF homologous factor receptor) (FGFR-like protein) (Fibroblast growth factor receptor 5) (FGFR-5), FGFRL1, FGFR5 FHFR
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FGFRL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.