



# GHRHR Antibody, HRP conjugated

Product Code	CSB-PA009413LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02643
Immunogen	Recombinant Human Growth hormone-releasing hormone receptor protein (23-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth hormone-releasing hormone receptor (GHRH receptor) (Growth hormone-releasing factor receptor) (GRF receptor) (GRFR), GHRHR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GHRHR