



GHRHR Antibody, Biotin conjugated

Product Code	CSB-PA009413LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02643
Immunogen	Recombinant Human Growth hormone-releasing hormone receptor protein (23-127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Growth hormone-releasing hormone receptor (GHRH receptor) (Growth hormone-releasing factor receptor) (GRF receptor) (GRFR), GHRHR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GHRHR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.