




lacZ Antibody

Product Code	CSB-PA009476LA01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00722
Immunogen	Recombinant Escherichia coli Beta-galactosidase protein (2-1024AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-galactosidase (Beta-gal) (EC 3.2.1.23) (Lactase), lacZ
Immunogen Species	Escherichia coli
Research Area	Tags & Cell Markers
Target Names	lacZ

Image



Western Blot

Positive WB detected in: Recombinant protein

All lanes: lacZ antibody at 7.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 117 kDa

Observed band size: 117 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.