



lacZ Antibody, HRP conjugated

Product Code	CSB-PA009476LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00722
Immunogen	Recombinant Escherichia coli Beta-galactosidase protein (2-1024AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-galactosidase (Beta-gal) (EC 3.2.1.23) (Lactase), lacZ
Immunogen Species	Escherichia coli (strain K12)
Research Area	Tags & Cell Markers
Target Names	lacZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.