



GRPR Antibody, FITC conjugated

Product Code	CSB-PA009943LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30550
Immunogen	Recombinant Human Gastrin-releasing peptide receptor protein (326-384AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gastrin-releasing peptide receptor (GRP-R) (GRP-preferring bombesin receptor), GRPR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRPR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.