



HAS1 Antibody, Biotin conjugated

Product Code	CSB-PA010139LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92839
Immunogen	Recombinant Human Hyaluronan synthase 1 protein (151-271AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hyaluronan synthase 1 (EC 2.4.1.212) (Hyaluronate synthase 1) (Hyaluronic acid synthase 1) (HA synthase 1) (HuHAS1), HAS1, HAS
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	HAS1