



HIP1R Antibody, HRP conjugated

Product Code	CSB-PA010366LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75146
Immunogen	Recombinant Human Huntingtin-interacting protein 1-related protein (304-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Huntingtin-interacting protein 1-related protein (HIP1-related protein) (Huntingtin-interacting protein 12) (HIP-12), HIP1R, HIP12 KIAA0655
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HIP1R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.