



HOXA4 Antibody, HRP conjugated

Product Code	CSB-PA010654LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q00056
Immunogen	Recombinant Human Homeobox protein Hox-A4 protein (23-173AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-A4 (Homeobox protein Hox-1.4) (Homeobox protein Hox-1D), HOXA4, HOX1D
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.