



HOXB4 Antibody, HRP conjugated

Product Code	CSB-PA010664LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17483
Immunogen	Recombinant Human Homeobox protein Hox-B4 protein (88-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-B4 (Homeobox protein Hox-2.6) (Homeobox protein Hox-2F), HOXB4, HOX2F
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXB4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.