



HOXD3 Antibody, FITC conjugated

Product Code	CSB-PA010686LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31249
Immunogen	Recombinant Human Homeobox protein Hox-D3 protein (263-346AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-D3 (Homeobox protein Hox-4A), HOXD3, HOX1D HOX4A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXD3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.