




HOXD8 Antibody

Product Code	CSB-PA010688LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13378
Immunogen	Recombinant Human Homeobox protein Hox-D8 protein (59-108AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-D8 (Homeobox protein Hox-4E) (Homeobox protein Hox-5.4), HOXD8, HOX4E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HOXD8

Image



Western Blot
 Positive WB detected in: PC-3 whole cell lysate,
 Mouse brain tissue
 All lanes: HOXD8 antibody at 4.5µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 32, 9 kDa
 Observed band size: 32 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.