



KCNK3 Antibody, HRP conjugated

Product Code	CSB-PA012071LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14649
Immunogen	Recombinant Human Potassium channel subfamily K member 3 protein (278-379AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel subfamily K member 3 (Acid-sensitive potassium channel protein TASK-1) (TWIK-related acid-sensitive K(+) channel 1) (Two pore potassium channel KT3.1) (Two pore K(+) channel KT3.1), KCNK3, TASK TASK1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNK3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.