



LAMA2 Antibody, Biotin conjugated

Product Code	CSB-PA012726LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24043
Immunogen	Recombinant Human Laminin subunit alpha-2 protein (974-1162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit alpha-2 (Laminin M chain) (Laminin-12 subunit alpha) (Laminin-2 subunit alpha) (Laminin-4 subunit alpha) (Merosin heavy chain), LAMA2, LAMM
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LAMA2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.