



LHX2 Antibody, Biotin conjugated

Product Code	CSB-PA012920LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50458
Immunogen	Recombinant Human LIM/homeobox protein Lhx2 protein (197-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM/homeobox protein Lhx2 (Homeobox protein LH-2) (LIM homeobox protein 2), LHX2, LH2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	LHX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.